

CD68 Polyclonal Antibody APC Conjugated

Catalog No: #C91466APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | CD68 Polyclonal Antibody APC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | Flow-Cyt ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 111-210/335 |
| Conjugates | APC |
| Target Name | CD68 |
| Other Names | CD68 |
| Accession No. | Uniprot:P31996Gene ID:12514 |
| Uniprot | P31996 |
| GeneID | 12514 |
| Calculated MW | 37 |
| Formulation | Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4 |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only