

CD68 Polyclonal Antibody FITC Conjugated

Catalog No: #C91467FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD68 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Ms Rat
Immunogen Type	151-250/335
Conjugates	FITC
Target Name	CD68
Other Names	CD68
Accession No.	Uniprot:P31996Gene ID:12514
Uniprot	P31996
GeneID	12514
Calculated MW	37
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only