

## VWF Polyclonal Antibody HRP Conjugated

Catalog No: #C91470HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	VWF Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	1701-1800/2813
Conjugates	HRP
Target Name	VWF
Other Names	VWF
Accession No.	Uniprot:P04275Gene ID:7450
Uniprot	P04275
GeneID	7450
Calculated MW	309
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only