

Semaphorin 5A Polyclonal Antibody FITC Conjugated

Catalog No: #C91472FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Semaphorin 5A Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	671-770/1074
Conjugates	FITC
Target Name	Semaphorin 5A
Other Names	Semaphorin 5A
Accession No.	Uniprot:Q13591Gene ID:9037
Uniprot	Q13591
GeneID	9037
Calculated MW	118
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only