

Sialoadhesin/CD169 Polyclonal Antibody PE Conjugated

Catalog No: #C91473PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sialoadhesin/CD169 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Ms Rat
Immunogen Type	1301-1400/1709
Conjugates	PE
Target Name	Sialoadhesin/CD169
Other Names	Sialoadhesin/CD169
Accession No.	Uniprot:Q62230Gene ID:20612
Uniprot	Q62230
GeneID	20612
Calculated MW	186
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only