

Cytokeratin 18 Polyclonal Antibody APC Conjugated

Catalog No: #C91474APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytokeratin 18 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	53-150/430
Conjugates	APC
Target Name	Cytokeratin 18
Other Names	Cytokeratin 18
Accession No.	Uniprot:P05783Gene ID:3875
Uniprot	P05783
GeneID	3875
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only