

Tau Polyclonal Antibody PE Conjugated

Catalog No: #C91475PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tau Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/758
Conjugates	PE
Target Name	Tau
Other Names	Tau
Accession No.	Uniprot:P10636Gene ID:4137
Uniprot	P10636
GeneID	4137
Calculated MW	52/79
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only