

LOX 1 Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C91476PE-Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LOX 1 Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-273/273
Conjugates	PE-Cy5
Target Name	LOX 1
Other Names	LOX 1
Accession No.	Uniprot:Q9XTA8Gene ID:100009322
Uniprot	Q9XTA8
GeneID	100009322
Calculated MW	31/50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only