

## BTN2A1 Polyclonal Antibody PE Conjugated

Catalog No: #C91477PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	BTN2A1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu
Conjugates	PE
Target Name	BTN2A1
Other Names	BTN2A1
Accession No.	Uniprot:Q7KYR7Gene ID:11120
Uniprot	Q7KYR7
GeneID	11120
Calculated MW	57
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only