

Leptin receptor Polyclonal Antibody FITC Conjugated

Catalog No: #C91479FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Leptin receptor Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms
Immunogen Type	561-660/1165
Conjugates	FITC
Target Name	Leptin receptor
Other Names	Leptin receptor
Accession No.	Uniprot:P48357Gene ID:3953
Uniprot	P48357
GeneID	3953
Calculated MW	103/132
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only