

iNos/Nos-2 Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C91485PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	iNos/Nos-2 Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1051-1153/1153
Conjugates	PE-Cy7
Target Name	iNos/Nos-2
Other Names	iNos/Nos-2
Accession No.	Uniprot:P35228Gene ID:4843
Uniprot	P35228
GeneID	4843
Calculated MW	130
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only