

NRXN3 Polyclonal Antibody PE Conjugated

Catalog No: #C91488PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NRXN3 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	28-120/1643
Conjugates	PE
Target Name	NRXN3
Other Names	NRXN3
Accession No.	Uniprot:Q9HDB5Gene ID:9369
Uniprot	Q9HDB5
GeneID	9369
Calculated MW	178
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only