

Caspase-1 P10 Polyclonal Antibody FITC Conjugated

Catalog No: #C91489FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-1 P10 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms
Immunogen Type	321-404/404
Conjugates	FITC
Target Name	Caspase-1 P10
Other Names	Caspase-1 P10
Accession No.	Uniprot:P29466Gene ID:834
Uniprot	P29466
GeneID	834
Calculated MW	10/45
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only