

N-cadherin Polyclonal Antibody Biotin Conjugated

Catalog No: #C91490Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-cadherin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	501-600/906
Conjugates	Bio
Target Name	N-cadherin
Other Names	N-cadherin
Accession No.	Uniprot:P19022Gene ID:1000
Uniprot	P19022
GeneID	1000
Calculated MW	100
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only