

OPG Polyclonal Antibody HRP Conjugated

Catalog No: #C91491HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	OPG Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/401
Conjugates	HRP
Target Name	OPG
Other Names	OPG
Accession No.	Uniprot:O00300Gene ID:4982
Uniprot	O00300
GeneID	4982
Calculated MW	44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only