

# Integrin beta 1/CD29 Polyclonal Antibody PE Conjugated

Catalog No: #C91493PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Integrin beta 1/CD29 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms
Immunogen Type	501-600/798
Conjugates	PE
Target Name	Integrin beta 1/CD29
Other Names	Integrin beta 1/CD29
Accession No.	Uniprot:P05556Gene ID:3688
Uniprot	P05556
GeneID	3688
Calculated MW	88
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only