

HMGB1 Polyclonal Antibody FITC Conjugated

Catalog No: #C91494FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HMGB1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	61-150/215
Conjugates	FITC
Target Name	HMGB1
Other Names	HMGB1
Accession No.	Uniprot:P09429Gene ID:100862258
Uniprot	P09429
GeneID	100862258
Calculated MW	25
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only