

Galectin 3 Polyclonal Antibody HRP Conjugated

Catalog No: #C91509HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galectin 3 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Rat
Immunogen Type	121-200/250
Conjugates	HRP
Target Name	Galectin 3
Other Names	Galectin 3
Accession No.	Uniprot:P17931Gene ID:3958
Uniprot	P17931
GeneID	3958
Calculated MW	29
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only