

5HT2A Receptor Polyclonal Antibody Biotin Conjugated

Catalog No: #C91517Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	5HT2A Receptor Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	171-270/471
Conjugates	Bio
Target Name	5HT2A Receptor
Other Names	5HT2A Receptor
Accession No.	Uniprot:P28223Gene ID:3356
Uniprot	P28223
GeneID	3356
Calculated MW	52
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only