

Nurr1 Polyclonal Antibody PE Conjugated

Catalog No: #C91518PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nurr1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	111-200/598
Conjugates	PE
Target Name	Nurr1
Other Names	Nurr1
Accession No.	Uniprot:P43354Gene ID:4929
Uniprot	P43354
GeneID	4929
Calculated MW	67
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only