

GPR49/LGR5 Polyclonal Antibody APC Conjugated

Catalog No: #C91519APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR49/LGR5 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	331-430/907
Conjugates	APC
Target Name	GPR49/LGR5
Other Names	GPR49/LGR5
Accession No.	Uniprot:O75473Gene ID:8549
Uniprot	O75473
GeneID	8549
Calculated MW	98
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only