

SERPIN A11 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91529Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SERPIN A11 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Ms Rat
Immunogen Type	181-180/427
Conjugates	Bio
Target Name	SERPIN A11
Other Names	SERPIN A11
Accession No.	Uniprot:Q8CIE0Gene ID:380780
Uniprot	Q8CIE0
GeneID	380780
Calculated MW	45
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only