

ompF Polyclonal Antibody HRP Conjugated

Catalog No: #C91533HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ompF Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC
Species Reactivity	Escherichia Coli Yersinia enterocolitica sub
Immunogen Type	171-213/379
Conjugates	HRP
Target Name	ompF
Other Names	ompF
Accession No.	Uniprot:A1JMN7Gene ID:4713258
Uniprot	A1JMN7
GeneID	4713258
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC=1:50-200

Note: This product is for in vitro research use only