

# Angiotensin II Type 1 Receptor Polyclonal Antibody PE Conjugated



Catalog No: #C91539PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Angiotensin II Type 1 Receptor Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/359
Conjugates	PE
Target Name	Angiotensin II Type 1 Receptor
Other Names	Angiotensin II Type 1 Receptor
Accession No.	Uniprot:P30556Gene ID:185
Uniprot	P30556
GeneID	185
Calculated MW	39
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only