

TEM1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91542Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TEM1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	251-350/757
Conjugates	Bio
Target Name	TEM1
Other Names	TEM1
Accession No.	Uniprot:Q9HCU0Gene ID:57124
Uniprot	Q9HCU0
GeneID	57124
Calculated MW	80
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only