

SDR-O Polyclonal Antibody Biotin Conjugated

Catalog No: #C91549Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | SDR-O Polyclonal Antibody Biotin Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IHC ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 181-280/313 |
| Conjugates | Bio |
| Target Name | SDR-O |
| Other Names | SDR-O |
| Accession No. | Uniprot:Q8NEX9Gene ID:121214 |
| Uniprot | Q8NEX9 |
| GeneID | 121214 |
| Calculated MW | 35 |
| Formulation | Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4 |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only