

RASSF10 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91552Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RASSF10 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/507
Conjugates	Bio
Target Name	RASSF10
Other Names	RASSF10
Accession No.	Uniprot:A6NK89 Gene ID:644943
Uniprot	A6NK89
GeneID	644943
Calculated MW	57
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only