

SLC15A3 Polyclonal Antibody FITC Conjugated

Catalog No: #C91554FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC15A3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	431-530/581
Conjugates	FITC
Target Name	SLC15A3
Other Names	SLC15A3
Accession No.	Uniprot:Q8IY34Gene ID:51296
Uniprot	Q8IY34
GeneID	51296
Calculated MW	64
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only