

Cyp2J3 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91556Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyp2J3 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Ms Rat
Immunogen Type	401-502/502
Conjugates	Bio
Target Name	Cyp2J3
Other Names	Cyp2J3
Accession No.	Uniprot:P51589Gene ID:313375
Uniprot	P51589
GeneID	313375
Calculated MW	55
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only