

## SLC36A2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91558FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SLC36A2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	261-360/483
Conjugates	FITC
Target Name	SLC36A2
Other Names	SLC36A2
Accession No.	Uniprot:Q495M3Gene ID:153201
Uniprot	Q495M3
GeneID	153201
Calculated MW	53
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only