

SLC39A13 Polyclonal Antibody HRP Conjugated

Catalog No: #C91560HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC39A13 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/371
Conjugates	HRP
Target Name	SLC39A13
Other Names	SLC39A13
Accession No.	Uniprot:Q96H72Gene ID:91252
Uniprot	Q96H72
GeneID	91252
Calculated MW	39
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only