

TNF alpha Polyclonal Antibody Cy5 Conjugated

Catalog No: #C91567Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF alpha Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu
Immunogen Type	181-235/233
Conjugates	Cy5
Target Name	TNF alpha
Other Names	TNF alpha
Accession No.	Uniprot:P01375Gene ID:7124
Uniprot	P01375
GeneID	7124
Calculated MW	17/26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only