

# Cytokeratin 19 Polyclonal Antibody PE-Cy7 Conjugated



Catalog No: #C91579PE-Cy7

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytokeratin 19 Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/403
Conjugates	PE-Cy7
Target Name	Cytokeratin 19
Other Names	Cytokeratin 19
Accession No.	Uniprot:P08727Gene ID:3880
Uniprot	P08727
GeneID	3880
Calculated MW	44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only