

## EGFP Polyclonal Antibody Biotin Conjugated

Catalog No: #C91580Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	EGFP Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Recombinant EGFP
Conjugates	Bio
Target Name	EGFP
Other Names	EGFP
Accession No.	Uniprot:P42212Gene ID:107331116
Uniprot	P42212
GeneID	107331116
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at -80°C

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only