

CD51 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91593Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD51 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/1048
Conjugates	Bio
Target Name	CD51
Other Names	CD51
Accession No.	Uniprot:P06756Gene ID:3685
Uniprot	P06756
GeneID	3685
Calculated MW	95/113
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only