

Mouse serum albumin Polyclonal Antibody Cy5 Conjugated

Catalog No: #C91594Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse serum albumin Polyclonal Antibody Cy5 Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | ICC IF |
| Species Reactivity | Hu Ms |
| Conjugates | Cy5 |
| Target Name | Mouse serum albumin |
| Other Names | Mouse serum albumin |
| Accession No. | Uniprot:P02769Gene ID:280717 |
| Uniprot | P02769 |
| GeneID | 280717 |
| Calculated MW | 68 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only