

H5N1 Hemagglutinin Polyclonal Antibody Biotin Conjugated

Catalog No: #C91597Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	H5N1 Hemagglutinin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Influenza A virus H5N1
Immunogen Type	200-300/568
Conjugates	Bio
Target Name	H5N1 Hemagglutinin
Other Names	H5N1 Hemagglutinin
Accession No.	Uniprot:O56140Gene ID:
Uniprot	O56140
Calculated MW	64
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only