

## CD162 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91600Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD162 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu
Immunogen Type	251-350/412
Conjugates	Bio
Target Name	CD162
Other Names	CD162
Accession No.	Uniprot:Q14242Gene ID:6404
Uniprot	Q14242
GeneID	6404
Calculated MW	41
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only