

## ZO-1 Polyclonal Antibody PE Conjugated

Catalog No: #C91608PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ZO-1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms
Immunogen Type	1451-1550/1748
Conjugates	PE
Target Name	ZO-1
Other Names	ZO-1
Accession No.	Uniprot:Q07157Gene ID:7082
Uniprot	Q07157
GeneID	7082
Calculated MW	191
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only