

## IBDV-VP2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91609Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IBDV-VP2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	IBDV
Immunogen Type	21-120/1012
Conjugates	Bio
Target Name	IBDV-VP2
Other Names	IBDV-VP2
Accession No.	Uniprot:P61825Gene ID:
Uniprot	P61825
Calculated MW	109
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only