

Ki-67 Polyclonal Antibody HRP Conjugated

Catalog No: #C91612HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ki-67 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Ms
Immunogen Type	1231-1330/2938
Conjugates	HRP
Target Name	Ki-67
Other Names	Ki-67
Accession No.	Uniprot:Q91VE6Gene ID:17345
Uniprot	Q91VE6
GeneID	17345
Calculated MW	358
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only