

ROR gamma T Polyclonal Antibody HRP Conjugated

Catalog No: #C91616HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | ROR gamma T Polyclonal Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | WB IHC ICC |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 1-100/516 |
| Conjugates | HRP |
| Target Name | ROR gamma T |
| Other Names | ROR gamma T |
| Accession No. | Uniprot:P51449Gene ID:19885 |
| Uniprot | P51449 |
| GeneID | 19885 |
| Calculated MW | 54 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:500-2000 IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only