

ROR gamma T Polyclonal Antibody APC Conjugated

Catalog No: #C91617APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ROR gamma T Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Ms Rat
Immunogen Type	1-100/516
Conjugates	APC
Target Name	ROR gamma T
Other Names	ROR gamma T
Accession No.	Uniprot:P51449Gene ID:19885
Uniprot	P51449
GeneID	19885
Calculated MW	54
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only