

CD163 Polyclonal Antibody APC Conjugated

Catalog No: #C91619APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD163 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/1156
Conjugates	APC
Target Name	CD163
Other Names	CD163
Accession No.	Uniprot:Q86VB7Gene ID:9332
Uniprot	Q86VB7
GeneID	9332
Calculated MW	121
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only