

Caveolin-1 Polyclonal Antibody HRP Conjugated

Catalog No: #C91621HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caveolin-1 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	2-80/178
Conjugates	HRP
Target Name	Caveolin-1
Other Names	Caveolin-1
Accession No.	Uniprot:Q03135Gene ID:857
Uniprot	Q03135
GeneID	857
Calculated MW	20
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only