

ITGA7 Polyclonal Antibody PE Conjugated

Catalog No: #C91626PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | ITGA7 Polyclonal Antibody PE Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 1081-1181/1181 |
| Conjugates | PE |
| Target Name | ITGA7 |
| Other Names | ITGA7 |
| Accession No. | Uniprot:Q13683Gene ID:3679 |
| Uniprot | Q13683 |
| GeneID | 3679 |
| Calculated MW | 25/125 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only