

ALK Polyclonal Antibody PE Conjugated

Catalog No: #C91628PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ALK Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	19-120/1620
Conjugates	PE
Target Name	ALK
Other Names	ALK
Accession No.	Uniprot:Q9UM73Gene ID:238
Uniprot	Q9UM73
GeneID	238
Calculated MW	174
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only