

CXCL10/IP10 Polyclonal Antibody HRP Conjugated

Catalog No: #C91631HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CXCL10/IP10 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC
Species Reactivity	Hu
Immunogen Type	29-70/98
Conjugates	HRP
Target Name	CXCL10/IP10
Other Names	CXCL10/IP10
Accession No.	Uniprot:P02778Gene ID:3627
Uniprot	P02778
GeneID	3627
Calculated MW	10
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only