

IGF2R/M6PR Polyclonal Antibody PE Conjugated

Catalog No: #C91632PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IGF2R/M6PR Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	701-800/2491
Conjugates	PE
Target Name	IGF2R/M6PR
Other Names	IGF2R/M6PR
Accession No.	Uniprot:P11717Gene ID:3482
Uniprot	P11717
GeneID	3482
Calculated MW	300
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only