

Hepatitis B Virus Surface Antigen Polyclonal Antibody Biotin Conjugated



Catalog No: #C91633Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatitis B Virus Surface Antigen Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hepatitis B virus(HBVD2, HBVC5)
Immunogen Type	31-100/389
Conjugates	Bio
Target Name	Hepatitis B Virus Surface Antigen
Other Names	Hepatitis B Virus Surface Antigen
Accession No.	Uniprot:P03142Gene ID:2703546
Uniprot	P03142
GeneID	2703546
Calculated MW	31/44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only